




## Prestan Professional Manikins

### Repairing a Faulty Infant Manikin Rate Monitor

**Step 1:** Using a flat head screwdriver, remove Monitor Cover as shown.




**Step 2:** Pull out Monitor.






**Step 3:** Remove batteries from Monitor and inspect each of the six wire contacts for corrosion. If mild corrosion is found, soak a Q-tip in isopropyl alcohol and cautiously brush the contacts to remove corrosion. If moderate corrosion is found, cautiously brush the contacts with a Scotch-Brite™ pad to remove corrosion. If severe corrosion is found, the monitor may need to be replaced.



**Step 4:** Inspect the two round gold contacts inside the manikin as shown for corrosion.





**Step 5A:** If mild corrosion or discoloration is found on either of the round gold contacts, soak a Q-tip in isopropyl alcohol and clamp in needle nose pliers. Cautiously brush the contacts to remove corrosion.



**Step 5B:** If moderate/severe corrosion or discoloration is found on either of the round gold contacts, cut a ½" strip from a Scotch-Brite™ pad and clamp in needle nose pliers. Cautiously brush the contacts to remove corrosion.



**Step 6:** Insert Monitor.



**Step 7:** Attach Monitor Cover.

