

**Adult Torso Skin
Replacement-
Medium Tone
RPP-ASKIN-MS**




**Adult Compression
Rate Feedback Monitor
RPP-AMON-1**




**Adult Head Assembly
RPP-AHEAD-1-MS**



**Ear Pin
Medium Tone
10138-MS**






Head Attachment Posts



Chest Frame Clip


**Adult Frame
RPP-AFRAME-1**

**Adult
Compression Rate
Feedback Monitor
RPP-AMON-1**




**Chest Plate
RPP-APLATE-
1**

**Adult Foam
Insert
RPP-AFOAM-1**




**Adult Frame
RPP-AFRAME-1**

Adult Frame Clip



**Adult Foam Insert
RPP-AFOAM-1**

Torso Skin Off



** MS = Medium Tone
DS = Dark Tone